Natrelle

UNDERSTANDING YOUR OPTIONS

Your consultation discussion guide

Selecting the right breast implant is an important part of the reconstruction process. Talk to your plastic surgeon about the look you want. Your surgeon can help you choose the best option based on your needs, wants, body structure, and more.



Natrelle® Breast Implants IMPORTANT SAFETY INFORMATION AND APPROVED USES

Breast implants are not considered lifetime devices. The longer people have them, the greater the chances are that they will develop complications, some of which will require more surgery.

Breast implants have been associated with the development of a cancer of the immune system called breast implant–associated anaplastic large cell lymphoma (BIA-ALCL). This cancer occurs more commonly in patients with textured breast implants than smooth implants, although rates are not well defined. Some patients have died from BIA-ALCL.

Patients receiving breast implants have reported a variety of systemic symptoms, such as joint pain, muscle aches, confusion, chronic fatigue, autoimmune diseases, and others. Individual patient risk for developing these symptoms has not been well established. Some patients report complete resolution of symptoms when the implants are removed without replacement.

QUESTIONS TO ASK YOUR PLASTIC SURGEON



Get to know your plastic surgeon	 How long have you been performi How many breast reconstructions 	astic surgery from a recognized and accredited program? ng breast reconstruction surgery? with implants do you perform a year? nd what is the most common reoperation you perform?
Understanding breast reconstruction surgery	using own body tissues—autologo 3. What are the risks and complicatio 4. What are the pros and cons for dire 5. Should I wait to have breast recon 6. What if my cancer recurs after bre 7. How long is breast reconstruction	ast reconstruction (implant-based reconstruction, bus, or implant combined with own body tissue)? ons of each type of breast reconstruction surgery? oct-to-implant and two-stage breast reconstruction surgery? struction until I have completed all my cancer treatments? ast reconstruction?
Understanding Natrelle® breast implants	 Do you have before-and-after brea what results are reasonable for me What do you take into account wh I want? 	s, will I need tissue expanders? What is the process? ast reconstruction surgery photos to help me understand a to expect? een choosing the best <i>Natrelle®</i> option to get the look lecide to choose <i>Natrelle®</i> breast implants?
What can I expect post surgery?	 What can I expect during my record. What are common risks and composition. What will my scars look like? How much pain or discomfort will When can I resume my normal act Do <i>Natrelle</i>[®] breast implants have 	lications associated with breast implants? I feel, and for how long? ivities?
Natrelle [®] Breast Implants IMPORTANT SAFETY INFORMATION AND APPROVED USES (continued)		 Breast reconstruction. This includes primary breast reconstruction to replace breast tissue that has been removed due to cancer or trauma or that has failed to develop properly due to a

Who can get breast implants?

Natrelle[®] Breast Implants are approved for the following:

- Breast augmentation for women at least 22 years old for silicone-filled implants and for women at least 18 years old for saline-filled implants. Breast augmentation includes primary breast augmentation to increase the breast size and revision surgery to correct or improve the result of a primary breast augmentation
- Breast reconstruction. This includes primary breast reconstruction to replace breast tissue that has been removed due to cancer or trauma or that has failed to develop properly due to a severe breast abnormality. This also includes revision surgery to correct or improve the result of a primary breast reconstruction



Natrelle[®] Breast Implants IMPORTANT SAFETY INFORMATION AND APPROVED USES (continued)

Who should NOT get breast implants?

Breast implant surgery should NOT be performed in:

- Women with active infection anywhere in their body
- Women with existing cancer or precancer of their breast who have not received adequate treatment for those conditions
- · Women who are currently pregnant or nursing

What should I tell my doctor?

Tell your doctor if you have any of the following conditions, as the risks of breast implant surgery may be higher:

- Autoimmune diseases (eg, lupus and scleroderma)
- A weakened immune system (eg, taking medications to decrease the body's immune response)
- Planned chemotherapy or radiation therapy following breast implant placement
- Conditions or medications that interfere with wound healing and blood clotting
- · Reduced blood supply to breast tissue
- Clinical diagnosis of depression or other mental health disorders, including body dysmorphic disorder and eating disorders
- Those with a diagnosis of depression or other mental health disorders should wait for resolution or stabilization of these conditions prior to undergoing breast implantation surgery

What else should I consider?

- There is a Boxed Warning for breast implants. Please see bold text at beginning
- Many changes to your breasts following implantation are irreversible. If you
 later choose to have your implants removed and not replaced, you may
 experience dimpling, puckering, wrinkling, or other cosmetic changes, which
 may be permanent
- Breast implantation is likely not a one-time surgery. The longer implants are in
 place, the greater the potential risk for complications. You will likely need additional
 surgeries on your breasts due to complications or unacceptable cosmetic results.
 Thus, you should also consider the complication rates for later (revision) surgery
 since you may experience these risks in the future

- Cancer treatments and surgery will affect the outcome and timing of breast reconstruction
- Breast implants may affect your ability to breastfeed, either by reducing or eliminating milk production
- Rupture of a silicone-filled breast implant is most often silent. Even if you
 have no symptoms, you should have your first ultrasound or MRI at 5 to 6 years
 after your initial implant surgery and then every 2 to 3 years thereafter regardless
 of whether your implants are for augmentation or reconstruction. If you have
 symptoms of or uncertain ultrasound results for breast implant rupture, an MRI
 is recommended. Additional imaging may be required depending on your medical
 history and status. The health consequences of a ruptured silicone gel-filled
 breast implant have not been fully established
- Routine screening mammography for breast cancer will be more difficult, and implants may rupture during the procedure. Perform self-examination every month for cancer screening and ask your surgeon to help you distinguish the implant from your breast tissue. Lumps, persistent pain, swelling, hardening, or changes in implant shape should be reported to your surgeon and possibly evaluated with imaging

What are key complications with breast implants?

Key complications include reoperation, implant removal with or without replacement, implant rupture with silicone-filled implants, implant deflation with saline-filled implants, and capsular contracture (severe scar tissue around the implant). Other complications include breast pain, swelling, asymmetry, wrinkling/rippling, implant malposition nipple complications, hypertrophic scarring, and implant palpability/visibility.

Talk to your doctor about other complications.

For more information, see the patient brochures at www.allergan.com/products.

To report a problem with *Natrelle*[®] Breast Implants, please call Allergan[®] at 1-800-624-4261.

The sale and distribution of *Natrelle*[®] Breast Implants is restricted to licensed physicians who provide information to patients about the risks and benefits of breast implant surgery.

Allergan Aesthetics an AbbVie company © 2022 AbbVie. All rights reserved. NATRELLE® and its design are trademarks of Allergan, Inc., an AbbVie company. Natrelle.com NAT64505-v4 05/22 014232

